

Biofeedback Driven



Biofeedback Driven™ is a new concept in biofeedback software. It gives you the opportunity to reduce monotony – hence increasing biofeedback training efficiency. BIOFEEDBACK DRIVEN™ animations give you a true choice in animation selection. Library I consists of 10 Cartoons & Puzzles, and 3 Slide Shows. You control which one is used and when. Biofeedback Driven™ is meant to be used in conjunction with Brainwave Animation Pro and 2.0 AVI unlock and enhances its motivational value. Each Cartoon holds the key to eager cooperation and learning – a child’s work is a child’s play. Each engaging full story-line cartoon embeds text and a moral. This series gives you a choice and it’s easy to move from cartoon to cartoon. When using a Cartoon you engage the child, as a willing participant, in the biofeedback process. Activating any animation is simple. Simply select the animation of your choice and begin the biofeedback session, as usual.

Each Puzzle requires the trainee remove the top layer of graphic to get to the animating reward. Younger kids particularly like this challenge. They are quick and simple, working well for children who have low tolerance for frustration and little ability to delay gratification.

The three Slide Shows transition from artistically rendered photo to photo, animating their way through each slide in the show. Adults and kids alike enjoy this series. **Here’s what your colleagues have said:** “The cartoons are great and the children love them. We have not encountered any technical problems. Thanks.” Carol Manchester, Ph.D. “I have used your animation a couple of times and am sold! Send me the set.” David Thomson, Ph.D. **Working with Kids or Adults, help them succeed!**

The animations are captivating – dramatically enhancing participation in treatment; The animations are bright and colorful; The Cartoons contain text, a moral and full story-line; The Puzzles entice the child to the second level where the animation lies; The Slide Shows intrigue the trainee with their mystifying transitions.

System Recommended/Required Processor: Pentium or like capability (586) Hard Drive: Recommend 1.5 GB available space (Less is workable) Memory: Recommend 32 Mb, or sufficient to run Windows on your PC. Biofeedback Software Recommended: BrainMaster systems using Brainwave Animation Pro or 2.0 AVI Viewer.

Bundled Price: \$575 (Affiliates \$475.00)

Dr. Padgitt’s Animation Libraries are shipped on CD, ready for easy installation. If you need technical assistance, Brainwave Center is there for you. The Brainwave Center is committed to an expanding series of animation Libraries to serve you and your clients.



BrainMaster Technologies, Inc.

24490 Broadway Avenue
 Oakwood Village, Ohio 44146
 Corporate Offices (440) 232-6000 Fax (440) 232-7171
 info@brainm.com